

CLOSURE PLUS PROTECTION



Clo-Sur P.A.D.™ and
Clo-Sur^{PLUS} P.A.D.™

Pressure Applied Dressings

Clo-Sur P.A.D.TM and Clo-Sur^{PLUS} P.A.D.TM

Following catheterization, achieving hemostasis with manual compression can be time-consuming and costly. Patients are also at increased risk for catheter-related blood stream infections (CRBSI)¹. The Clo-Sur P.A.D. accelerates hemostasis. The Clo-Sur^{PLUS} P.A.D. antimicrobial dressing protects skin surface puncture sites against nosocomial infection.



Clo-Sur^{PLUS} P.A.D. provides 6 day antimicrobial barrier protection on skin surface puncture sites.




CLOSURE PLUS PROTECTION

CLOSURE

Non-woven topical pad made with a proprietary polypropylate, a positively-charged material with known hemostatic activity that expedites clot formation.

Faster hemostasis leads to decreased hold time and time to ambulation.

Can also be used for patients on anticoagulation therapy.



¹ <http://www.safecarecampaign.org/CRBSI.html>

² *Eur Radiol*, 2007 Mar 17(3):693-700.

³ WuXiAppTec Report 867961 May 31, 2011.

Unlike other vascular closure devices, the Clo-Sur P.A.D. and Clo-Sur^{PLUS} P.A.D. have no contraindications and can be used by patients with known sensitivity to bovine-derived materials.

Water soluble for painless and easy removal, which is more comfortable for the patient.

PLUS PROTECTION

Clo-Sur^{PLUS} P.A.D. provides an antimicrobial barrier for up to 6 days.

Can be used as antimicrobial wrap around indwelling lines to protect skin surface puncture sites.

MICROBIAL ASSAYS TEST RESULTS³

ORGANISM	7 DAY INCUBATION		
	REPLICATE 1	REPLICATE 2	REPLICATE 3
A. brasiliensis	No Growth	No Growth	No Growth
B. cereus	No Growth	No Growth	No Growth
C. albicans	No Growth	No Growth	No Growth
E. coli	No Growth	No Growth	No Growth
E. faecium	No Growth	No Growth	No Growth
P. aeruginosa	No Growth	No Growth	No Growth
S. aureus	No Growth	No Growth	No Growth
S. pyogenes	No Growth	No Growth	No Growth

Distributed by:



Catalog Number	Product Description
MVP44	4x4 cm Hemostasis Pad
SCV22P-AB	4x4 cm Hemostasis Pad with Antimicrobial Barrier
HPAD22	Kit with SCV22P-AB and a 3M TM Tegaderm TM Film*

* 3M and Tegaderm are trademarks of 3M
Packaged 10 per box.

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT
 Merit Medical Europe, Middle East, and Africa (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 22
 Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733
 Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30
 Finland 0800 770 586 • Germany 0800 182 0871 • Ireland (Republic) 1800 553 163 • Italy 800 897 005 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84
 Norway 800 11629 • Sweden 020 792 445 • UK 0800 973 115

www.merit.com